

All-new!

Pixel-perfect acuity testing.



**ClearChart® 4, 4X, and 4P**

Intuitive, connected, and purpose-built in the USA.

**Reichert**  
TECHNOLOGIES  
**AMETEK**

Advancing Eye Care. American Innovation.



## ClearChart® 4 Family of Digital Acuity Systems

Pixel-perfect acuity testing · [reichert.com/clearchart](http://reichert.com/clearchart)

The all-new family of ClearChart® Digital Acuity Systems from Reichert® are all-in-one, elegant, simple-to-use, and uniquely designed to meet all of your acuity testing needs with a wide selection of optotypes and other special test charts. The **ClearChart 4 Digital Acuity System**, **ClearChart 4X Enhanced Digital Acuity System**, and **ClearChart 4P Polarized Digital Acuity System** all feature a new 24 inch, high resolution, LED backlit display, are lightweight, cool running and long-lasting. ClearChart also seamlessly connects to a Reichert digital Phoropter® for the best possible refraction experience for both you and your patient.

	ClearChart 4	ClearChart 4X	ClearChart 4P
Purpose-built for acuity testing	•	•	•
Wide range of optotypes and tests	•	•	•
24" high resolution LED backlit display	•	•	•
Cool running, free of moving fans or hard drives	•	•	•
Pediatric fixation cartoon video	•	•	•
Connects to Reichert digital Phoroptors	•	•	•
Astigmatic testing	•	•	•
ETDRS		•	•
Contrast sensitivity		•	•
Patient education slides		•	•
Import your own videos and images		•	•
Stereo testing			•
Color vision testing			•



**AMETEK®** © 2017 AMETEK, Inc. & Reichert, Inc. (3-2017) · [www.reichert.com](http://www.reichert.com) ·   

ClearChart® is made in the USA. Reichert, Reichert Technologies, Phoropter, and ClearChart are registered trademarks of Reichert, Inc.